

# Full Cast High Noble Crown & Bridge

**CDT Code: D2790**



## Features

- Long-term wear characteristics
- Full cast crowns can last for 30-40 years with good oral hygiene
- Full cast gold is still one of the most wear friendly of dental materials and is extremely bio-compatible
- Superior biocompatibility: gold crowns have been shown to have fewer pathogens due to highly polished, thin margin areas and the fact that the pathogens do not have an affinity for gold

## Benefits

- Longevity
- Plaque resistance
- High marginal integrity and durability

## Cementation

- Cementation using standard cementation materials & technique

## Contraindications

- Not indicated for anterior restorations where high esthetics are required

## Indications

- Indicated for single unit crowns, multi-unit bridges, inlays and onlays

## Prep

- Occlusal/Incisal Reduction: 1.5mm
- Facial/Lingual Reduction: 1.2 mm
- Cervical Reduction: 1.0 mm reduction; however, a knife edge marginal prep is acceptable



Tampa, Florida 33634  
877-337-7800 | www.DDSLlab.com

\*\*The person signing this authorization and/or the dental practice accepts responsibility for payment of the related charges and agrees to pay all legal and collections costs in the event the account is in collections or litigation, including reasonable fees. LAB-300670 © 2015 DDS Lab. All rights reserved. DL10334, 3274, L005, 552

## Full Cast High Noble RX - Laboratory Procedure Authorization

### REQUIRED INFORMATION

This is my first case with DDS Lab

Doctor Name \_\_\_\_\_  
Last First

Practice Name \_\_\_\_\_

Address \_\_\_\_\_

Phone \_\_\_\_\_

E-mail \_\_\_\_\_

Patient Name \_\_\_\_\_

Patient Chart # \_\_\_\_\_  M  F Age \_\_\_\_\_

Rx Date \_\_\_\_\_ Due Date/Delivery by 5pm on \_\_\_\_\_

**(standard working time if no date given)**

*Case turnaround times are based on the date the prescription is received at DDS Lab. Please allow at least 10 business days (M-F) from that date.*

### CASE INSTRUCTIONS

**Tooth Numbers** Please circle single units and bracket splinted units.

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
32	31	30	29	28	27	26	25	24	23	22	21	20	19	18	17

#### Metal

Full Cast HN (AU ≥40)

#### Restoration

Crown  
 Bridge

#### Return for

Die trim  Finish

### CROWN DESIGN

#### Pontic Design



#### If Insufficient Room

Trim opposing\*  Reduction coping

#### Occlusal Clearance

Light\*  Open  Tight

\*Standard design if an option is not selected.

#### Contact

Light\*  Medium  Heavy

### SPECIFIC INSTRUCTIONS

Please provide any photos, study models, diagnostic casts with case.

\_\_\_\_\_

Dentist signature\*\*

Dentist license no.

Script has been reviewed for accuracy, legibility and completion.  
Impressions have been approved by the dentist.

### LAB USE ONLY

Impression \_\_\_\_\_ Model \_\_\_\_\_

Bite \_\_\_\_\_ Articulator \_\_\_\_\_ Parts \_\_\_\_\_

Crown  Bridge  Other

Appliance Enclosed \_\_\_\_\_