

Digital impressions - we work with all!



3Shape Trios

- Use the "add connections" feature to send request to cadcam@ddslab.com.
- We will get the request within 5 minutes. Once approved the doctor will be able to send scans.



Allied Star

- Log in to AS Connect Web.
- Click "My Partner" menu.
- Click the Add Partner button from the My Partner menu.
- Search for DDS LAB via Zip Code: "33634" or email address: cadcam@ddslab.com.
- Select DDS Lab and click "Apply".



Carestream

- Send request to cadcam@ddslab.com through "Settings" or "Configuration".
- Files can be sent through Carestream portal or upload the .STL file through DDS portal.



DDX

- DDX user should search for DDS Lab in Florida
- Verify the phone number associated: 877-337-7800.



Dentsply Sirona

- Use "add favorite lab" feature to add DDS Lab to your favorites.
- Search for DDS LAB via Zip Code: "33634".
- Scan can be sent right away- simply choose DDS Lab from the drop-down menu.



Dexis:

- Navigate to "Settings" or "Configuration".
- Send request using cadcam@ddslab.com.



E4D

- DDS Lab receives E4D scans through DDX (Henry Schein portal).
- Doctor can email the STL file or they can upload through DDS or DDX portals.



iTero

- Search for "DDS Lab" or by Company ID (9866) through "add preferred lab". If either option is unavailable, please call iTero directly to add DDS Lab as a preferred lab.
- When the doctor scans any case, the file goes to Align Tech and iTero team will evaluate the scan. If the scan needs to be taken again, iTero will contact the doctor. If the scan is approved, it will be sent to the lab within 8 hours.
- Ortho files must be uploaded through iRecord in the iTero portal. Files will still be accepted through the iTero case managers but that will delay the case.



Ivoclar:

- Navigate to "Settings" or "Configuration".
- Send request using cadcam@ddslab.com.



Medit Link

- Login to your Medit Link account
- Use "search partners" to find DDS Lab by name or email address: cadcam@ddslab.com.



NeoScan:

- Navigate to "Settings" or "Configuration".
- Send request using cadcam@ddslab.com.



Shining3D

- Use "Connections" to access "+NEW" tab.
- Search "DDS Lab" via Zip Code: 33634.
- Ensure to select account with "TOOTH" image.