

Lab-Fabricated. Clinically Guided. Designed for Real-World Removable Workflows.

Standardized for consistency. Flexible for clinical need.

Our Platform Approach

TrueDent™

Default for most indications.

Lucitone Digital Print™

Used selectively based on clinical need or doctor preference.

Why We Offer Both

Different clinical situations require different 3D printed denture solutions.



TrueDent™ (Default)

Used for:

- Immediate dentures
- Copy / replacement dentures
- High-volume, repeatable workflows

Why:

- Monolithic (no tooth-to-base interface)
- Consistent outcomes
- Efficient and scalable

Lucitone Digital Print™

Used for:

- Implant conversion or overdenture workflows
- Cases requiring alignment with existing dentures or partials
- Multi-component prosthetic cases

Why:

- Supports coordinated prosthetic workflows
- Enables cross-arch material and tooth matching
- Adds flexibility for specific clinical needs