

## Digital Workflow Information

### For Photos:

- Before submitting a case/scan via Dexis IS, make a note on the Rx notepad that photos will be sent separately for combination with the case.
- Submit the scan as usual.
- Once the scan has been submitted, email the photos to [DDSLABPIX@DDSLAB.COM](mailto:DDSLABPIX@DDSLAB.COM) to ensure they are properly attached to the case.

### For Digital Remakes:

Shade Change Remakes:

- Log in to MYDDSLAB Portal and create an online Rx.
- Download the Rx and email it to [DigitalRemake@ddslab.com](mailto:DigitalRemake@ddslab.com).
- We will reuse the original design file and refabricate the case in the updated shade.

### Margin Issue Remakes:

- We recommend rescanning the prep to confirm the margin issue was not caused by lack of clarity in the original scan.
- On the Rx notepad, note that the case is a remake and provide a brief summary of the issue.

### Angel Aligners STL File Export:

Currently, Angel Aligners integration is not available. However, you can export STL files by following these steps:

1. Open the scan in DEXIS IS ScanFlow.
2. Select the One-click Export option.
3. Choose the folder where you want to save the file.
4. Select STL as the file type